



OneBAT Training on Bat Sanitary Surveillance

20-21 May 2026 | Padova (Italy)

**The Development of Pseudotyped Viruses for Cell Tropism
and Serosurveillance Studies of Bat Viruses**

Simon D Scott

s.d.scott@kent.ac.uk

University of
Kent




"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

Bat Virus Detection



- The identification of virus infection in animals (including bats) and man is usually via detecting viral genomic material (DNA, RNA) or proteins
- A positive test result tells us that the bat either has a current infection or is recovering from it
- The 'window of opportunity' to detect these molecules is relatively narrow – days or weeks usually
- For wild animals this would mean sampling would have to take place during this period
- An alternative detection method is to identify **virus-specific antibodies**, which have been produced during the immune response to the infection
- Antibody production begins during the infection and gradually builds
- These antibodies circulate for months or years ('memory')
- This gives a wider 'window of opportunity' for detection




Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

Antibody (Serology) Tests



- There are a number of methods of detecting specific antibodies – serological tests
- Some measure the **binding** of antibodies to virus-specific proteins, often those on the surface of the virus, which are the target of immunological responses in the host
- Others measure anti-virus biological activity – antibody **neutralisation** of viruses
- Neutralisation assays can be conducted using native, replicative virus
- However if investigating pathogenic or potentially pathogenic viruses, these assays have to be carried out in high biological containment facilities for safety reasons
- Such facilities are only found in a relatively few institutes/laboratories, which can act as a 'bottleneck' for research
- However, there is an alternative approach to perform pathogenic virus antibody neutralisation assays – using so-called **pseudotyped viruses** as surrogates





Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

Pseudotyped Viruses (PVs) & Utilisation

- Pseudotyped viruses (PVs) are hybrid particles with a protein **'core' (capsid)** encasing a virus genetic material (genome), encapsulated in an **'envelope'** derived from host cell outer membrane
- 'Pseudotyping' occasionally occurs naturally after co-infection of cells by 2 different virus strains, releasing offspring (progeny) with mixed components derived from the 2 parents
- Importantly, the **envelope** displays surface proteins from the other (heterologous) virus
- These surface proteins are **essential for virus infection** (attachment & entry) and **stimulate host immune responses (e.g. antibodies)**
- PVs can be artificially engineered, and are usually non-replicative – thus need low biosafety containment
- PVs can be utilised to **study virus entry, cell tropism, receptors, immune responses (antibody neutralisation) and antiviral screening**





PVs have 3 main components:


- **Envelope** glycoprotein(s)
- **Core** (e.g. derived from a lentivirus)
- **Transgene** (e.g. 'reporter gene' for quantifying infection or inhibition of infection (i.e. immunity))



Pseudotyped Virus Generation



- PVs can be generated by the introduction of essential genes/components on small circular DNA molecules called **plasmids** (via genetic engineering) into a **'producer' cell line**
 - Once the plasmids enter the cells; expression of genes into proteins, synthesis of new genomes then assembly of virus **'core'** particles within cell
 - Finally, the PV particles escape from the producer cells by **'budding'** (pushing through outer cell membrane) and being wrapped in the host-derived **envelope** (containing specific virus glycoproteins - **GPs**)
 - The success of the process (*which is by no mean guaranteed!*) is measured via a **'titration' assay** – determines concentration of biologically active viral particles
 - If concentration or **'titre'** is sufficient the PVs can be used for antibody detection
- **Pseudotyped Virus Neutralisation Tests/Assays (PVNT/PVNA)**




Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

Pseudotyped Virus (PV) Generation & Titration

3 PV components:



Use of PVs to investigate cell 'tropism'



- PVs can be utilised to examine which cells can be infected by study viruses – i.e. cells from different animals and humans (some mammal lines developed in OneBAT)
- This enables a wide range of cells to be analysed without need to handle the native pathogen
- For OneBAT we investigated the 'tropism' of 3 different bat viruses in a number of cell types: Lloviu filovirus (LLOV) & 2 lyssaviruses - **Lleida & West Caucasian Bat Viruses (LLEBV & WCBV)**
- All 3 viruses could infect the hamster, dog and human cells we tested, with varying efficiency
- Interestingly, only cells from certain bat species could be infected
- Included cells from *Miniopterus schreibersii* (known host of the 3)
- Most permissive cells were human kidney (HEK293)
- Thus HEK293 utilised as target cells in future assays





Funded by the
European Union


"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

Pseudotyped Virus Neutralisation Test/Assay (PVNT/PVNA)

- Without presence of antibodies PVs can enter target cells and display the matching **'receptor'** on their surface → **luminescence**
- If PVs are incubated with a serum sample containing antibodies, these antibodies interfere with the PV entering target cells → luminescence reduced ↓
- Thus if a bat has had a prior infection with a virus, these antibodies will 'neutralise' the virus → luminescence reduced ↓
- Amount of luminescence inversely proportional to the amount of antibodies (VNABs) in bat serum



Screening bat serum samples for viruses



- The **Viral Pseudotype Unit (VPU)**, University of Kent, UK has over a decade of experience in generating and utilising many different PVs for research
- VPU collaborations with some members of the OneBAT consortium before grant initiated – i.e. working on the **Lloviu filovirus (LLOV)** with Dr Kemenesi @ University of Pécs, Hungary
- As part of the OneBAT project ~1400 bats (from multiple roosts in 5 countries, over several years) have had small volumes of blood taken and sent to VPU
- Screening of sera via PVNA for antibodies to **LLOV, LLEBV & WCBV** – thus each sample analysed for VNABs from 3 different viruses
- For this **‘Triplex’** system we need 3 separate, independent **reporter** read-outs
- Thus, we generated 3 different PVs & 3 separate assays – then combined



PV Neutralisation Assays – reporter genes


- PVNA carried out in 96 well assay plates
- Serial dilution of bat sera + standard PV amount
- Initial approach was using Green Fluorescent protein (GFP) gene in LLOV PVs, Firefly Luciferase (Fluc) for WCBV & Renilla Luciferase (Rluc) for LLEBV
- These reporters can all be detected independently without interference – 3 different readouts for the 3 virus/antibody combinations

1. Add standard sera mix and perform 2-fold dilutions

2. Add pseudotypes mix

3 Target cells

LLOV +
LLEBV +
WCBV PVs



Aequorea victoria



Lampyridae spp.



Renilla reniformis




Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."


PVNA - data readouts



- GFP read via microscopy or image analyser (clear plate)
- In VPU done by manually counting PV infected, green cells
- Less green cells if antibodies present



4. Read the plate fluorescence



5. Add Dual-Glo® Luciferase Assay Reagent to the plate.



20-25 °C 10 minutes

6. Transfer into a white plate and measure firefly luminescence



7. Add Dual-Glo® Stop & Glo® Reagent to the plate.



20-25 °C 10 minutes

8. Measure *Renilla* luminescence



Luciferase read using a Microplate Reader
- Fluc then Rluc (different substrates to produce luminescence)



"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."

OneBAT serum sample 'triple' screening – ongoing...



- Green (GFP) cell counting proved too time consuming manually
- Would need expensive image analyser to automate
- For practical reasons now testing LLOV separately (using Fluc), with lyssaviruses still Fluc & Rluc
- Readout appear as a **'heat map'** with numerical values to analyse
- Experiments include positive & negative control sera
- **Screening of the ~1400 bat sera ~50% complete**



Funded by the
European Union

“Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them.”

Serum screening data analysis – ongoing...



- Transfer data to spreadsheets & analysis of all the data being produced
- Comparison of test sample data to controls
- Statistical analysis
- Compiling study data to gain overall picture of virus distribution across Europe

Controls:

Cell only

PV only

Negative serum

Positive serum

Triplex WCBV FLuc s601-624 Quick Read 2026.02.27 19:53:01

WCBV

Luminescence None, 0.3Sec

All	1	2	3	4	5	6	7	8	9	10	11	12
A	50	60	10	70	1.83E+04	1.05E+04	5.21E+04	9.53E+04	80	80	1.6E+05	1.22E+05
B	130	160	190	170	1.55E+05	1.23E+05	1.81E+05	2.19E+05	1.2E+04	4.64E+03	2.78E+05	1.72E+05
C	2.42E+05	3.51E+05	1.37E+05	1.49E+05	140	120	150	210	1.35E+05	1.08E+05	170	70
D	3.25E+05	3.77E+05	2.93E+05	2.64E+05	1.08E+04	40	3.18E+04	5.73E+04	2.58E+05	2.46E+05	1.16E+04	2.44E+04
E	2.65E+05	3.05E+05	4.37E+03	2.12E+03	100	40	50	170	9.04E+04	8.75E+04	4.4E+04	7.92E+04
F	3.22E+05	4.29E+05	2.25E+05	1.41E+05	4.3E+03	210	460	140	2.57E+05	2.45E+05	3.04E+05	2.18E+05
G	210	250	1.87E+03	90	40	30	30	50	130	170	590	8.25E+03
H	4.99E+04	3.82E+04	6.94E+04	4.11E+04	1.34E+04	8.5E+03	5.35E+03	530	1.11E+04	7.06E+03	1.2E+05	1.13E+05

Triplex LLEBV RLuc Dual-Glo Luciferase Assay System 2026.02.27 20:18:04

LLEBV

Luminescence - Renilla luciferase activity measurement None, 0.3Sec

All	1	2	3	4	5	6	7	8	9	10	11	12
A	80	170	5.06E+05	4.28E+05	140	70	80	1.02E+03	1.75E+05	1.5E+05	3.97E+05	2.14E+05
B	270	360	6.6E+05	7.81E+05	180	90	1.18E+05	1.97E+05	7.64E+05	6.45E+05	7.71E+05	5.43E+05
C	6.99E+05	1.09E+06	4.76E+04	3.25E+04	4.75E+03	1.63E+04	1.95E+05	2.1E+05	3.77E+05	5.19E+05	2.82E+05	3.03E+05
D	1.03E+06	1.03E+06	4.66E+05	4.25E+05	3.89E+05	1.45E+05	5.35E+05	8.65E+05	7.84E+05	5.8E+05	8.61E+05	7.26E+05
E	9.62E+05	6.99E+05	1.25E+03	190	170	260	1.58E+05	1.74E+05	3.2E+05	4.09E+05	4.52E+04	1.07E+05
F	9.85E+05	1.42E+06	4.61E+05	2.64E+05	5.23E+03	430	6.37E+05	6.56E+05	9.67E+05	9.12E+05	5.95E+05	5.01E+05
G	5E+05	2.97E+05	310	5.2E+03	8.64E+05	7.79E+05	3.41E+05	4.44E+05	430	350	1E+05	1.07E+05
H	4.85E+05	6.27E+05	9.06E+04	8.24E+04	1.03E+06	9.13E+05	8.84E+05	7.59E+05	8.31E+04	5.7E+04	5.28E+05	3.98E+05



Funded by the
European Union

“Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them.”



Simon Scott

University of Kent

s.d.scott@kent.ac.uk



Maria Francesca Priore

IZSVE

MFPriore@izsvenezie.it



Edward Wright

University of Sussex

e.wright@sussex.ac.uk



**Thanks to all OneBAT
samplers and cell line
developers!**



Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."



Thanks for your attention!



www.onebat.eu



Funded by the
European Union

"Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Health and Digital Executive Agency (HADEA). Neither the European Union nor the granting authority can be held responsible for them."